

# ***ClearBand* 10x Tris-Glycine Buffer**

500 ml, 1000 ml

**Cat No:** TG05, TG10

**Shipping** : Ship at ambient temperature.

**Storage** : Store at 4-8°C.

## **General Information**

***ClearBand*** Tris-Glycine Buffer (10x) is an ideal stock solution for preparing standard Tris-glycine transfer buffer used for Western blotting. ***ClearBand*** Tris-Glycine Buffer (10x) is prepared with ultra pure water and filter sterilized.