

***ClearBand* 10% SDS**

100 ml, 500 ml

Cat No: S100, S500

Shipping : Ship at ambient temperature.

Storage : Store at room temperature.

General Information

ClearBand 10% SDS (10% w/v) is sodium dodecyl sulfate prepared in 18 MΩ water. SDS is a detergent that is used for denaturation of proteins especially in polyacrylamide gel electrophoresis applications.

Precautions

If precipitation forms, incubate the solution at 50°C with occasional mixing until precipitate is completely dissolved.

Do not store *ClearBand* 10% SDS solution below 15°C.