

## ***ClearBand* Potassium Acetate**

**3M, pH: 5.5**

**Ready-to-use Solution**

100 ml

**Cat No:** PA100

**Shipping** : Ship at ambient temperature.

**Storage** : Store between 2-8°C.

### **General Information**

*ClearBand* Potassium Acetate, 3M pH:5.5 Ready-to-use Solution is prepared using molecular biology grade potassium acetate in ultrapure water. *ClearBand* Potassium Acetate, 3M pH:5.5 Ready-to-use Solution is 0.22 µm filter sterilized and suitable for purification and precipitation of nucleic acids. It can be used for also precipitating dodecyl sulfate (DS) and DS-bound proteins allowing the removal of proteins from DNA.