

EcoZyme RNase A

1 ml, 5x1 ml

Cat No: EZRA1, EZRA5

Shipping : Ship at ambient temperature.

Storage : Store at -20°C.

General Information

EcoZyme RNase A is a pre-made molecular biology grade enzyme solution which is supplied at a concentration of 10 mg/ml in 10 mM Tris-HCl (pH 8). *EcoZyme* RNase A from bovine pancreas is a stable and highly reactive, DNase and protease-free, endoribonuclease that specifically degrades single stranded RNA at C and U residues. *EcoZyme* RNase A is pretreated to remove DNase I and is suitable for RNA digestion in plasmid and genomic DNA purification procedures.