

***ClearBand* 10x Tris-Buffered Saline with Tween[®] 20, pH:7.6**

500 ml, 1000 ml

Cat No: TBST05, TBST10

Shipping : Ship at ambient temperature.

Storage : Store at room temperature. Use within 1 week once it is diluted to 1x.

General Information

ClearBand 10x Tris-Buffered Saline with Tween[®] 20 (TBST) is a balanced pH stabilizing salt solution used for especially western blotting and ELISA procedures. *ClearBand* TBS enables washing without disturbing antibody-antigen binding interactions.

ClearBand 10x TBST should be diluted to 1x working solution. 1x formulation contains 0.1% Tween 20.

ClearBand 10x TBST is prepared with ultra pure water and filter sterilized.