

***ClearPeak* Oligo(dT)₁₈**

50 rxn, 100 rxn

Cat No: Q1850, Q18100

Shipping : Ship with blue ice.

Storage : Store at -20°C. Avoid freeze and thaw cycles.

General Information

ClearPeak Oligo(dT)₁₈ primer is a string of 18 deoxythymidylic acid residues that hybridizes to the poly(A) tail of most eukaryotic mRNA and can be used as a primer for first strand cDNA synthesis with reverse transcriptase. The primer is supplied in distilled water. *ClearPeak* Oligo(dT)₁₈ primer is recommended for use with *EcoScript* 5x First Strand cDNA Synthesis Kit (Cat. no. ERT50, ERT100, ERT200).

We recommend using 1 µl of *ClearPeak* Oligo(dT)₁₈ primer per 20 µl cDNA synthesis reaction volume.

Related Products

ClearPeak 2x SYBR Master Mix (SMM1, SMM5)

EcoScript 5x First Strand cDNA Synthesis Kit (ERT50, ERT100, ERT200)