

***ClearBand* 10x Phosphate Buffered Saline, pH:7.4**

500 ml, 1000 ml, 5x500 ml

Cat No: PBS10x-500, PBS10x-1000, PBS10x-5005

Shipping : Ship at ambient temperature.

Storage : Store at room temperature.

General Information

ClearBand 10x Phosphate-Buffered Saline (PBS) is a balanced salt solution which is used for a number of cell culture applications including washing cells before trypsinization, transport of cells or tissue samples, diluting cells for counting, and preparing reagents. *ClearBand* PBS does not contain calcium and magnesium for rinsing chelators from the culture before cell dissociation.

ClearBand 10x PBS is prepared with ultra-pure water, filter sterilized and autoclaved.