

# **NutriCulture Antibiotic Antimycotic Solution (100×)**

100 ml

**Cat No:** AA100

**Shipping** : Ship with dry ice.

**Storage** : Store at -20°C.

## **General Information**

**NutriCulture** Antibiotic Antimycotic Solution (100x) is an antibiotic mixture that can be effectively used to prevent gram-positive and gram-negative bacterial and fungal contamination in cell culture media. **NutriCulture** Antibiotic Antimycotic Solution is provided as a ready to use 100x concentrate and recommended final concentration is 10 ml/L.

## **Components**

Penicillin G Sodium ( $10^7$  Units/L); Streptomycin Sulfate (10000 mg/L); Amphotericin B (25 µg/ml); NaCl (9000 mg/L).