

## ***NutriCulture* PBS 1x**

Phosphate Buffered Saline, w/o Magnesium Chloride, w/o Calcium Chloride  
Cell Culture Reagent

500 ml, 10x500 ml

**Cat No:** PBS1x-500, PBS1x-50010

**Shipping** : Ship at ambient temperature.

**Storage** : Store at room temperature.

### **General Information**

*NutriCulture* Phosphate-Buffered Saline (PBS) is a balanced salt solution, which is used for a number of cell culture applications including washing cells before trypsinization, transport of cells or tissue samples, diluting cells for counting, and preparing reagents. *NutriCulture* PBS does not contain calcium and magnesium for rinsing chelators from the culture before cell dissociation. *NutriCulture* PBS is prepared with ultra-pure water, filter sterilized and autoclaved.