

***NutriCulture* Penicillin Streptomycin Solution (100x)**

100 ml

Cat No: PS100

Shipping : Ship with dry ice.

Storage : Store at -20°C.

General Information

NutriCulture Penicillin Streptomycin Solution (100x) is an antibiotic mixture that can be effectively used to prevent Bacterial contamination in cell culture media. This mixture can be used to prevent contaminations with gram-positive and gram-negative bacteria. *NutriCulture* Penicillin Streptomycin Solution is provided as a ready to use 100x concentrate and recommended final concentration is 10 ml/L.

Components

Penicillin G Sodium (10^7 Units/L); Streptomycin Sulfate (10000 mg/L); NaCl (9000 mg/L)