

# **NutriCulture Trypsin-EDTA (0.25%)**

100 ml

**Cat No:** TE025

**Shipping** : Ship with dry ice.

**Storage** : Store at -20°C.

## **General Information**

**NutriCulture** Trypsin-EDTA (0.25%) contains EDTA and phenol red. **NutriCulture** Trypsin-EDTA (0.25%) is made from trypsin powder, which is a mixture of proteases derived from porcine pancreas. Trypsin-EDTA is widely used for cell culture applications including cell dissociation, routine cell culture passaging, and primary tissue dissociation.